



Tech Tracks™

Application of Alternative Methods of Body Temp Measurement in Swine

Vol III

Purpose: Evaluate use of LifeChip® as a means to track body temperature in swine


Investigators: Matt Allerson, Bob Morrison
College of Veterinary Medicine, University of Minnesota, St Paul MN

Test Methods

- Thirty pigs injected subcutaneously and intramuscularly with LifeChip BioThermo® RFID transponder
- Pigs housed at University of Minnesota isolation facility
- Three groups of 10 pigs
- A pig infected with H1N1 influenza was introduced into each group day 1
- Body temp was measured for seven consecutive days by four methods
 - Sub-Q transponder (SQ)
 - Intramuscular transponder (IM)
 - Infrared camera (IR)
 - Rectal thermometer (RT)
- Mean temps measurements were compared and correlated
- Mean temps validity measurements were calculated

Test Results

- Measurement Results
 - RT highest at 104.3 F
 - IM next at 103.9 F
 - SQ at 100.8 F
 - IR at 100.6 F
- RT and IM measurements were highly correlated
 - $r=0.86$, 95% bootstrap C.I. 0.79, 0.90
- IR and IM measurements were moderately correlated
 - $r=0.61$, 95% bootstrap C.I. 0.47, 0.69



Discussion

- Elevated body temp is useful predictor of disease
- Obtaining body temp measurements rectally can be time consuming and stressful
- Alternative methods like transponders or IR cameras can provide rapid results
- IM transponders were highly correlated to core body temp
- IM transponders seem ideal for monitoring body temp in research setting
- Further research is needed for review in other species and production settings



Leading Since 1945
www.destronfearing.com

Destron Fearing™

© 2011 Destron Fearing
950-0079-001